

**Care and  
Cleaning guide  
Upholstery**

**Vinterno**



## Disinfectant use:

The information published in this care guide refers to the performance of Vinterno nv upholstery of products in specifics tests conducted under laboratory conditions. Results may vary under actual conditions.

- Asepticare TB + II
- Antifect (in 0.5 % dilution)
- Caviwipes XL Packed
- Cleanisept Wipes
- Clorox Healthcare Bleach Germicidal Wipes
- Coverage Plus Germicidal Wipes
- Desomed rapid FD 366
- Hydrogen Peroxide Cleaner Disinfectant Wipes
- Instrunet Anionsurf
- Lysol Foaming Disinfectant Cleaner
- Microzid
- Oxivir TB Wipes
- Oxycide Diluted -3oz/gl in water
- PDI Sani-cloth bleach wipes
- Sanicloth AF 3
- Sodium Hypochlorite. Disinfectant Bletach cleaner. Dilution: 1:10 mix of bleach and water.
- Surfanios Premium (in 0,25% dilution)
- Surfa'Safe Premium (in 0,25% dilution)
- Super Sani-cloth
- Virex II 256 (1:256 Dilution)
- Virox 5 Rtu
- Wet wipe
- Wex Cide 128 (1:128 dilution)

Vinterno nv is not sponsored by, affiliated with, or endorsed by the above trademark owners. Additionally, any use by Vinterno nv of any third-party trademark does not constitute any endorsement by Vinterno nv of such third parties.

Rinse the area with clean water, then dry with a clean, lint-free cloth after using any disinfectants and cleaners.

## Cleaning after disinfection:

- The vinyl should be cleaned periodically in order to maintain its appearance and prevent build-up of dirt and contaminants.
- Any stain, spills, or soiling should be cleaned up immediately to prevent the possibility of permanent staining.
- Use a mixture of liquid soap (neutral pH) and water or our approved cleaning products for coated fabrics to remove stains on the surface of the material. Remove only with a damp white cloth.
- Lacquers, strong cleaners, detergents, xylene-based solutions, acetone, or ketone (MEK) cause immediate damage and contribute to the deterioration of the material. The use of such cleaners is at the owner's risk.

## Daily cleaning:

1. Clean the entire surface using a mixture of 1:9 liquid soap (neutral pH) and water. Do not use soap or cleaning solutions that contain alcohol, ketones, xylene, acetates or solvents (mineral/white spirits).
2. Completely remove excess cleaner with a white, clean, damp cloth.
3. Dry the surface.

## Cleaning methods to remove tougher stains:

1. **Coffee, juice, red wine, sauces, chocolate, grease, crayon and suntan lotion.**  
Remove excess spill with a damp cloth. Clean with a 1:9 mix of neutral liquid soap and water. Next, rinse with clean water and dry.
2. **Make-up**  
Use a straight application of neutral PH soap. Wipe with a clean white cloth or soft brush, rinse with clean water, and dry.
3. **Blood, urine, feces**  
Use a 1:1 mix of Isopropyl alcohol and water. Rinse with clean water and dry.
4. **All other tough stains**  
Use a 0.5% solution of household bleach and water (1:9 mix of 5% concentrated household bleach and water). Rinse thoroughly with clean water and pat surface dry.

Note: For cleaning use a soft white cotton cloth saturated with the recommended cleaning agent and rub the stain in circles 10 times. Pat dry with another soft white cotton cloth and check the results. Do not use steel bristle brushes, abrasive scouring pads or metal, these cause permanent damage to the fabric.

Methods 2, 3 & 4 are not recommended for everyday cleaning!

After cleaning, rinse the area with clean water, then dry with a clean, lint-free cloth. Failure to remove chemical residue will cause discoloration of the coated fabric surface to manifest earlier than usual in the life cycle of the product.



Vinterno nv

Vermeirenplein 8  
2920 Kalmthout  
Belgium

+32 (0)3 640 32 80

[www.vinterno.com](http://www.vinterno.com)  
[info@vinterno.com](mailto:info@vinterno.com)